

R20 – Docking Arm Matrix



XS = 206 mm
S = 236 mm
L = 276 mm

Sitzhöhe hinten (mm)

	380	390	400	410	420	430	440	450	460	470	480	490	500	510	520	540	560
320						XS	XS	XS	XS	S	S	S	S	S	S	S	S
330						XS	XS	XS	XS	S	S	S	S	S	S	S	S
340					XS	XS	XS	XS	S	S	S	S	S	S	S	S	S
350					XS	XS	XS	XS	S	S	S	S	S	S	S	S	S
360				XS	XS	XS	XS	XS	S	S	S	S	S	S	S	S	S
370				XS	XS	XS	XS	XS	S	S	S	S	S	S	S	S	S
380			XS	XS	XS	XS	XS	S	S	S	S	S	S	S	S	S	S
390			XS	XS	XS	XS	XS	S	S	S	S	S	S	S	S	S	S
400			XS	XS	XS	XS	S	S	S	S	S	S	S	S	S	S	L
410			XS	XS	XS	XS	S	S	S	S	S	S	S	S	S	S	L
420						S	S	S	S	S	S	S	S	S	S	L	L
430						S	S	S	S	S	S	S	S	S	S	L	L
440				S	S	S	S	S	S	L	L	L	L	L	L	L	L
450				S	S	S	S	S	S	L	L	L	L	L	L	L	L
460			S	S	S	S	S	S	L	L	L	L	L	L	L	L	L
470			S	S	S	S	S	S	L	L	L	L	L	L	L	L	L
480		S	S	S	S	S	S	L	L	L	L	L	L	L	L	L	L
490		S	S	S	S	S	S	L	L	L	L	L	L	L	L	L	L
500		S	S	S	S	S	L	L	L	L	L	L	L	L	L	L	L
510		S	S	S	S	S	L	L	L	L	L	L	L	L	L	L	L
520		S	S	L	L	L	L	L	L	L	L	L	L	L	L	L	L
530		S	S	L	L	L	L	L	L	L	L	L	L	L	L	L	L
540		S	S	L	L	L	L	L	L	L	L	L	L	L	L	L	L
550		S	S	L	L	L	L	L	L	L	L	L	L	L	L	L	L
560		S	S	L	L	L	L	L	L	L	L	L	L	L	L	L	L
570		S	S	L	L	L	L	L	L	L	L	L	L	L	L	L	
580		S	S	L	L	L	L	L	L	L	L	L	L	L	L	L	
590		S	S	L	L	L	L	L	L	L	L	L	L	L	L	L	
600		L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	