


# R20 – Docking Arm Matrix



XS = 206 mm
S = 236 mm
L = 276 mm

		Sitzbreite (mm)																																
		A	B	320	330	340	350	360	370	380	390	400	410	420	430	440	450	460	470	480	490	500	510	520	530	540	550	560	570	580	590	600		
A: Hintere Sitzhöhe (mm)	380	115																																
	390	125																			S	S	S	S	S	S	S	S	S	S	S	S	L	
	400	135									XS	XS	XS	XS					S	S	S	S	S	S	S	S	S	S	S	S	S	S	L	
	410	145						XS	XS	XS	XS	XS	XS				S	S	S	S	S	S	S	S	L	L	L	L	L	L	L	L	L	
	420	155				XS	XS	XS	XS	XS	XS	XS	XS				S	S	S	S	S	S	S	L	L	L	L	L	L	L	L	L	L	
	430	165	XS	XS	XS	XS	XS	XS	XS	XS	XS	XS	XS		S	S	S	S	S	S	S	S	S	L	L	L	L	L	L	L	L	L	L	
	440	175	XS	XS	XS	XS	XS	XS	XS	XS	S	S	S	S	S	S	S	S	S	S	S	S	L	L	L	L	L	L	L	L	L	L	L	L
	450	185	XS	XS	XS	XS	XS	XS	S	S	S	S	S	S	S	S	S	S	S	S	L	L	L	L	L	L	L	L	L	L	L	L	L	L
	460	195	XS	XS	S	S	S	S	S	S	S	S	S	S	S	S	S	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L
	470	205	S	S	S	S	S	S	S	S	S	S	S	S	S	S	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L
	480	215	S	S	S	S	S	S	S	S	S	S	S	S	S	S	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L
	490	225	S	S	S	S	S	S	S	S	S	S	S	S	S	S	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L
	500	235	S	S	S	S	S	S	S	S	S	S	S	S	S	S	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L
	510	245	S	S	S	S	S	S	S	S	S	S	S	S	S	S	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L
	520	255	S	S	S	S	S	S	S	S	S	S	S	S	S	S	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L
	540	275	S	S	S	S	S	S	S	S	S	S	S	S	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L
560	295	S	S	S	S	S	S	S	S	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	L	